


# Extraction Forceps UK/European Patterns

## Extraction Forceps UK / European Patterns

Nova Extraction Forceps are manufactured from surgical grade stainless steel with a hinged joint to ensure optimum distribution of forces along the instrument. They are available in most common America, English and European patterns. If the forceps that you require are not shown we will be pleased to provide you with a quotation upon request.



| Fig. | Application                           | 65      | Upper premolars                        | 101/101A | Upper roots and bicuspid            |
|------|---------------------------------------|---------|--|----------|-------------------------------------|
| 1    | Upper centrals and canines            | 67      | Upper wisdoms                          | 107      | Upper canines                       |
| 2    | Upper laterals                        | 73      | Lower molars                           | 108      | Upper roots                         |
| 7    | Upper anteriors, premolars and roots  | 74/ 74M | Lower roots                            | 113      | Upper roots                         |
| 10   | Upper molars                          | 74N     | Small lower roots and crowded incisors | 123      | Children's lower teeth              |
| 13   | Lower bicuspid                        | 75      | Lower bicuspid                         | 136      | Children's upper centrals and roots |
| 16   | Lower molars                          | 75A     | Lower centrals and canines             | 137      | Lower centrals and roots            |
| 17   | Upper molars, right                   | 76      | Upper roots                            | 138      | Children's upper teeth and roots    |
| 18   | Upper molars, left                    | 76N     | Small upper roots                      | 143      | Lower roots                         |
| 19   | Upper wisdoms                         | 79      | Lower wisdom                           | 147      | Upper roots                         |
| 22   | Lower molars                          | 86      | Lower molars (Cowhorn)                 | 150      | Upper incisors                      |
| 23   | Lower molars                          | 87      | Lower molars (Cowhorn)                 | 151      | Lower incisors                      |
| 24   | Lower molars                          | 88      | Upper wisdom                           | 157      | Children's upper teeth              |
| 29   | Upper roots                           | 89      | Upper wisdom                           | 158      | Children's upper teeth              |
| 30   | Upper roots                           | 90      | Upper molars, right                    | 159      | Children's upper teeth              |
| 33   | Lower roots                           | 92      | Upper molars, left                     | 160      | Children's lower teeth              |
| 37   | Children's upper incisors and canines | 93      | Upper molars, right                    | 161      | Children's lower teeth              |
| 39   | Children's upper molars               | 94      | Upper molars, left                     | 162      | Children's lower teeth              |
| 44N  | Upper roots                           | 95      | Upper molars, right                    | 163      | Children's upper teeth              |
| 51   | Upper roots                           | 96      | Upper molars, left                     | 150      | Upper incisors                      |
| 52   | Upper roots                           | 97      | Small upper molars, right              | 151      | Lower incisors                      |
| 53   | Upper molars                          |         | Small upper molars, left               |          |                                     |

# Extraction Forceps UK/European Patterns



Extraction Forcep 1  
N0941



Extraction Forcep 2  
N0943



Extraction Forcep 7  
N0945



Extraction Forcep 17  
N0947



Extraction Forcep 18  
N0949



Extraction Forcep 19  
N0950



Extraction Forcep 22  
N0951



Extraction Forcep 23  
N0952



Extraction Forcep 24  
N0942



Extraction Forcep 29  
N0953



Extraction Forcep 30  
N0946



Extraction Forcep 33A  
N0948



Extraction Forcep 44N  
N1044



Extraction Forcep 51  
N0954




Extraction Forcep 51A  
N0956



Extraction Forcep 52  
N0958




Extraction Forcep 67  
N0959



Extraction Forcep 73  
N0961



Extraction Forcep 73S  
N0960




Extraction Forcep 74  
N0963



Extraction Forcep 74M  
N0962



Extraction Forcep 74N  
N0965



Extraction Forcep 75  
N0966



Extraction Forcep 75A  
N0964



Extraction Forcep 76  
N0967

# Extraction Forceps UK/European Patterns




Extraction Forcep 76N  
N0969



Extraction Forcep 76S  
N0968



Extraction Forcep 79  
N1920



Extraction Forcep 86  
N0971



Extraction Forcep 87  
N0970



Extraction Forcep 89  
N0978



Extraction Forcep 90  
N0980



Extraction Forcep 92  
N0982



Extraction Forcep 93  
N0984



Extraction Forcep 94  
N0972



Extraction Forcep 95  
N0974



Extraction Forcep 96  
N0986



Extraction Forcep 97  
N0988



Extraction Forcep 101  
N0990



Extraction Forcep 101A  
N0992



Extraction Forcep 107  
N0994




Extraction Forcep 108  
N0996



Extraction Forcep 113  
N0998



Extraction Forcep 136  
N1000



Extraction Forcep 137  
N1002



Extraction Forcep 143  
N1004



Extraction Forcep 147  
N1006


# Extraction Forceps UK/European Patterns



Eagle Beak Lower Univ.  
N1018



Eagle Beak Upper Left  
N1020



Eagle Beak Upper Right  
N1022




Forcep Rongeurs 3  
N1050



Forcep Rongeurs 4  
N1051



Rongeur Friedmann  
N1704



Rongeur Mini Friedmann  
N0757-Straight



Rongeur Mini Friedmann  
N0759-Curved



# Extraction Forceps UK/Childs & Pedodontic



Extraction Forcep 37  
N0955 Childs



Extraction Forcep 39  
N0957 Childs



Extraction Forcep 123  
N0973 Childs



Extraction Forcep 138  
N0975 Childs



Extraction Forcep 157  
N1008 Childs



Extraction Forcep 158  
N1010 Childs



Extraction Forcep 159  
N1012 Childs



Extraction Forcep 160  
N0976 Childs



Extraction Forcep 161  
N0977 Childs



Extraction Forcep 162  
N1014 Childs




Extraction Forcep 163  
N1016 Childs



Extraction Forcep 29S  
N0944 Childs and small  
upper roots



Pedodontic Forcep 1  
N1055




Pedodontic Forcep 2  
N1056



Pedodontic Forcep 3  
N1057




Pedodontic Forcep 4  
N1058



Pedodontic Forcep 5  
N1059



Pedodontic Forcep 6  
N1060



Pedodontic Forcep 7  
N1061